

## AMBULATORY AIDS



### Adult Walkers

- Anodized 1"aluminum frame with reinforced steel cross brace
- Simple 1 button release mechanism
- Foldable for easy transport or storage
- Soft, comfortable hypalon handgrips
- Anti-rattle silencers
- Available in 3 different sizes



### Aluminum Canes

- Anodized aluminum shaft
- Round, offset or orthopedic handles
- Soft, comfortable hypalon handgrips
- Adjustable in 1" increments



### Quad Canes

- Anodized aluminum shaft with reinforced steel base
- Soft, comfortable hypalon handgrips
- Handle rotates for left or right hand use
- Adjustable in 1" increments

## AMBULATORY AIDS

Encore Medical offers an assortment of affordable ambulatory aids by Tuffcare. Please see individual specification information below.

### AMBULATORY PRODUCT SPECIFICATIONS:

#### Walkers

Three sizes: Junior, Adult and Heavy Duty

Weight Limit: 300 lbs.

Note: Bariatric model available 400 lb capacity

#### Aluminum Canes

29" 40" High

Weight  $\frac{3}{4}$  lb

Weight Limit: 260 lbs.

#### Quad Canes

28" – 37" High

Weight  $1\frac{3}{4}$  lb

Weight Limit: 250 lbs.



12804 S. Spring Street  
Los Angeles, CA 90068  
(800)\_800-558-8636 Fax (310) 538-9177