

P R O D U C T S P E C I F I C A T I O N S

BORG WARNER ALPHA III



Featuring:

- Patient Hand Control
- Nurse Control Pedals
- $\frac{3}{4}$ Length Side Rails



BORG WARNER ALPHA III

The Borg Warner Alpha III bed is a great fit in today's Long-Term Care setting. Basic functionality and sturdy construction provides comfort and ease of use for both patient and caregiver alike. The heavy duty Nylasar plastic side rails provide complete safety and can collapse with a single lever function. This bed includes new Avalon style head & footboards and the Encore Medical 1 Year Warranty.

Includes: Folding $\frac{3}{4}$ Length Side rails, Footend Pedal Staff Controls with Patient Lockout, Pendant Hand Control, 5" Brake Casters, Pan Deck surface and New Head & Footboards.

Options: IV Pole, Full Length or Split Side rails, Wall Bumpers, 4" Casters

BORG WARNER ALPHA III DIMENSIONS:	BORG WARNER ALPHA III FUNCTIONS:
Length – 88"	Head Up / Down
Width (Rails up) – 39.5"	Hi-Low Elevation
Width (Rails down) – 37"	Knee Up / Down
Highest – 30"	Trendelenberg & Reverse Trendelenberg
Lowest – 15"	
Sleep Surface – 35 x 80"	
Weight Capacity - 500 lbs.	



12804 S. Spring Street
Los Angeles, CA 90068
(800) 800-558-8636 Fax (310) 538-9177