

P R O D U C T      S P E C I F I C A T I O N S

---

## HILL ROM 835/837



*Featuring:*

- Patient Rail Controls
- Tuck-away Side Rails
- Central Brake & Steer System

## HILL ROM 835/837

The Hill Rom 835/837 bed is our most popular bed. The functions include head, knee and hi-low elevation along with cardiac chair and Trendelenberg positions. The four split side rails provide complete safety and can tuck away underneath the frame with a single push function. The 835/837 can be controlled with dual side Central Brake & Steer Pedals and offers Automatic Contour at the knee section. This bed includes new Avalon style head & footboards and the Encore Medical 1 Year Warranty.

**Includes:** Central Brake & Steer, Automatic Contour, Split Side rails with Controls, Foot End Nurse Lockout Controls, 5" Casters, Spring or Pan Deck Surface and New Head & Footboards.

**Options:** IV Pole, Full Length Side rails, ¾ Length Side rails, Wall Bumpers, Trapeze Bracket, Pendant Hand Control, 3" or 4" Casters, Side rail Communication

HILL ROM 835/837 DIMENSIONS:	HILL ROM 835/837 FUNCTIONS:
Length – 88"	Head Up / Down
Width (Rails up) – 41.5"	Hi-Low Elevation
Width (Rails down) – 36.5"	Knee Up / Down
Highest – 30"	Cardiac Chair
Lowest – 18"	Trendelenberg & Reverse Trendelenberg
Sleep Surface – 35 x 80"	
Weight Capacity - 500 lbs.	
Under bed Clearance – 7.5"	



12804 S. Spring Street  
Los Angeles, CA 90068  
(800) 800-558-8636 Fax (310) 538-9177