

P R O D U C T S P E C I F I C A T I O N S

HILL ROM 840/842



Featuring:

- Patient Rail Controls
- Tuck-away Side Rails
- Central Brake & Steer System
- Retractable Frame

HILL ROM 840/842

The Hill Rom 840/842 bed is a great bed for either Long-Term or Acute Care. The functions include head, knee and hi-low elevation along with cardiac chair and Trendelenberg positions. The four split side rails provide complete safety and can tuck away underneath the frame with a single push function. The 840/842 bed also has a retractable frame that keeps the head section in line with the bedside cabinet for better patient comfort. This bed includes new Avalon style head & footboards and the Encore Medical 1 Year Warranty.

Includes: Retractable Frame, Central Brake & Steer, Automatic Contour, Split side rails with Controls, Foot End Nurse Lockout Controls, 5" Casters, Spring or Pan Deck Surface and New Head & Footboards.

Options: IV Pole, Wall Bumpers, Plastic Head & Footboards, Trapeze Bracket, Pendant Hand Control, Side rail Communication.

HILL ROM 840/842 DIMENSIONS:	HILL ROM 840/842 FUNCTIONS:
Length – 91"	Head Up / Down
Width (Rails up) – 41.5"	Hi-Low Elevation
Width (Rails down) – 36.5"	Knee Up / Down
Highest – 31"	Cardiac Chair
Lowest – 17"	Trendelenberg & Reverse Trendelenberg
Sleep Surface – 35 x 84"	Retractable Frame
Weight Capacity - 500 lbs.	
Under bed Clearance – 7.5"	

