

P R O D U C T S P E C I F I C A T I O N S

HILL ROM 850/852



Featuring:

- Patient Rail Controls
- Tuck-away Side Rails
- Central Brake & Steer System
- Retractable Frame
- Patient Nightlight
- CPR Release

HILL ROM 850/852

The Hill Rom 850/852 bed is an ideal fit in an Acute Care facility. The functions include head, knee and hi-low elevation along with cardiac chair and Trendelenberg positions. The four split side rails provide complete safety and can tuck away underneath the frame with a single push function. The 850/852 bed has a retractable frame that keeps the head section in line with the bedside cabinet for better patient comfort. Added benefits include a Patient Nightlight and CPR release. This bed includes new Avalon style head & footboards and the Encore Medical 1 Year Warranty.

Includes: Patient Nightlight, CPR Release, Retractable Frame, Central Brake & Steer, Automatic Contour, Split Side rails with Controls, Foot End Nurse Lockout Controls, 5" Casters, Spring or Pan Deck Surface and New Head & Footboards.

Options: IV Pole, Wall Bumpers, Plastic Head & Footboards, Trapeze Bracket, Pendant Hand Control, Side rail Communication.

| HILL ROM 850/852 DIMENSIONS: | HILL ROM 850/852 FUNCTIONS: |
|------------------------------|---------------------------------------|
| Length – 91" | Head Up / Down |
| Width (Rails up) – 41.5" | Hi-Low Elevation |
| Width (Rails down) – 36.5" | Knee Up / Down |
| Highest – 30.5" | Cardiac Chair |
| Lowest – 16.5" | Trendelenberg & Reverse Trendelenberg |
| Sleep Surface – 35 x 84" | Retractable Frame |
| Weight Capacity - 500 lbs. | |
| Under Bed Clearance – 7.5" | |



12804 S. Spring Street
Los Angeles, CA 90068
(800)_558-8636 Fax (310) 538-9177