

P R O D U C T S P E C I F I C A T I O N S

HILL ROM AFFINITY SERIES BIRTHING BED

Models Include:

Affinity, Affinity II & Affinity III



Featuring:

- Central Brake & Steer
- One Step Removable Foot Rest
- Tuckaway Side rails
- Adjustable Built-in Foot Supports



HILL ROM AFFINITY II

The Hill Rom Affinity Series birthing bed has been the industry workhorse in birthing units throughout the world. The functions include head, knee and hi-low elevations. This bed includes a refinished wood headboard and the Encore Medical 1 Year Warranty.

Includes: Placenta Basin, Calf Supports, Central Brake & Steer System, Nightlight, Lockout Controls, Instant CPR Release and New Custom fit Foam Mattress.

Options: IV Pole, Wall Bumpers, Side rail Communication, Pendant Hand Control, Air Mattress.

HILL ROM AFFINITY II DIMENSIONS:	HILL ROM AFFINITY II FUNCTIONS:
Length – 90”	Head Up / Down
Width (Rails up) – 42”	Hi-Low Elevation
Width (Rails down) – 34”	Knee Up / Down
Highest – 38”	Cardiac Chair
Lowest – 22”	Trendelenberg & Reverse Trendelenberg
Sleep Surface – 34 x 78”	
Weight Capacity - 500 lbs.	



12804 S. Spring Street
Los Angeles, CA 90068
(800) 800-558-8636 Fax (310) 538-9177