

P R O D U C T S P E C I F I C A T I O N S

HILL ROM CRITICAL CARE 894/896



Featuring:

- Patient Rail Controls
- Tuck-away Side rails
- Central Brake & Steer System
- CPR Release
- Radiolucent Head & Foot Section
- Patient Scale
- Electric or Manual Operation

HILL ROM CRITICAL CARE

The Hill Rom Critical Care bed is ideal for Intensive Care or Critical Care units. The functions include head, knee and hi-low elevation along with cardiac chair and hand controlled Trendelenberg positions. The four split side rails provide complete safety and can tuck away underneath the frame with a single push function. Radiolucent head section allows for in-bed X-Ray and this bed also includes a patient scale. The Critical Care bed also includes a CPR release, plastic head & footboards and the Encore Medical 1 Year Warranty.

Includes: Patient Scale, Radiolucent Head & Foot Section, Patient Nightlight, CPR Release, Central Brake & Steer, Split Side rails with Controls, Foot End Nurse Controls, 8" Casters, Pan Deck Surface and Plastic Head & Footboards.

Options: Air Mattress, IV Pole, Wall Bumpers, Pendant Hand Control, Side rail Communication.

HILL ROM CRITICAL CARE DIMENSIONS:	HILL ROM CRITICAL CARE FUNCTIONS:
Length – 91"	Head Up / Down
Width (Rails up) – 38.5"	Hi-Low Elevation
Width (Rails down) – 34"	Knee Up / Down
Highest – 40.5"	Cardiac Chair
Lowest – 22.5"	Trendelenberg & Reverse Trendelenberg
Sleep Surface – 33 x 80"	Electric & Manual Operation
Weight Capacity - 500 lbs.	
Under Bed Clearance – 7.5"	



12804 S. Spring Street
Los Angeles, CA 90068
(800) 800-558-8636 Fax (310) 538-9177