

MATTRESSES

CARE MATTRESS



Featuring:

- Soft Heel Pressure Reducing Section
- Dense Foam Side Rails Support
- Moisture Resistant Nylon Cover



CARE MATTRESS

The Encore Medical Care Mattress is an ideal foam mattress with a unique soft heel section for pressure reduction. The dense foam side rails provide additional support for moving in and out of bed. The cover is made of an Anti-Microbial nylon taffeta with rubber butyl backing for outstanding moisture resistance. Please see individual specification information below. This product carries a 5 Year Warranty.

CARE MATTRESS SPECIFICATIONS:

Cover System: Top Cover – Moisture Resistant, Anti-Microbial 70 denier green nylon with rubber back. Bottom Cover – Slide Resistant Anti-Microbial vinyl laminate with handles, French closure and two vents.

Foam Core: three layers of anti-microbial high-density foam with very firm foam side rails that add additional support for ingress/egress. 250 lb. capacity.

Optional: CTB 603 Fire-Barrier available. (California Required)

Standard sizes:

35" x 80"

35" x 84"

Custom Sizes Available



12804 S. Spring Street
Los Angeles, CA 90068
(800) 800-558-8636 Fax (310) 538-9177