

MATTRESSES

MULTIFLEX REVERSIBLE MATTRESS



Featuring:

- Reversible Foam Mattress
- Dense Foam Side Rails
- Moisture Resistant Nylon Cover
- 7 Year Warranty

MULTIFLEX MATTRESS

The Encore Medical Multiflex Mattress is a reversible foam mattress allowing it to be flipped to either side for support. The dense foam side rails provide additional support moving in and out of bed. The cover is made of an Anti-Microbial nylon taffeta with rubber butyl backing for outstanding moisture resistance. The foam surface is a flat cut convoluted foam for patient comfort and durability. Please see individual specification information below. This product carries a 7 Year Warranty.

MULTIFLEX SPECIFICATIONS:

Cover System: Moisture Resistant, Anti-Microbial 70 denier green nylon taffeta with rubber butyl backing. French closure and two vents. No handles on this model.

Foam Core: Three layers including anti-microbial foam with very firm foam side rails for additional support. The top and bottom layers of the mattress are 1.5" thick 1.85 lb. high density polyurethane foam, flat cut convoluted for patient comfort and durability. The center layer is a very firm 3" thick 1.9 lb. density for very firm support. Firm side rails through all layers provide extra support for ingress/egress. 250 lb capacity.

Optional: CTB 603 Fire-Barrier available. (California Required)

Standard sizes:

35" x 80"

35" x 84"

Custom Sizes Available



12804 S. Spring Street
Los Angeles, CA 90068
(800)_800-558-8636 Fax (310) 538-9177