




MAXXUM SERIES - LOW BEDS


MAXXUM Series - Electric



- Hand control
- Individual motor lockout
- Foot end controls



- 1/2 head, 1/2 foot or 3/4 length side rails
- Assist bar
- Side rail pads




- Perforated pan or orthopedic grid deck



- 4 IV and 4 traction receptacles



QD 8100 series - 77/8" - 28"



- MX lock system



- Total lock caster system




- High angle vascular frame

Electric Low Beds

30" high
 36" height
 19" height
 19" height
 19" height
 Affordable price



• Main frame mounted trapeze



• Bed extender and offset head/foot board brackets
 • Bed extender pad




• Bed travels virtually straight up and down



Foot rest



Transport/storage dollies



• DC motors and control
 • Battery back-up and charger

MAXXUM Series - Electric Low Beds



QD 8250 series - **Fast MAXXUM** - 7 7/8" - 30"

Specifications:

QD 8100 series - QD 8250 series

• Length including bumper	88"	88"
• Length (deck)	80"	80"
• Width (deck)	35"	35"
• Deck height - low	7 7/8"	7 7/8"
Deck height - high	28"	30"
Weight capacity	500 lbs.	500 lbs.

• *NOTE: Shorter beds also available*



A wide variety of pressure reduction mattresses for MAXXUM series beds.

Specifications and dimensions are approximate and subject to change without notice. Patents pending.

9/04

Encore Medical Technologies, Inc.

12804 S Spring Street, Los Angeles, CA 90061
 (800) 558-8636, (310) 538-9138, (310) 538-9177 fax
www.encoremedical.com, email: info@encoremedical.com

Introducing the NEW



MAXXUM Series - Fast Rising Beds



QD8250 with 360° mobility option

Features

- Fast - 7 7/8" low to 30" high in 20 seconds
- Mobile or locked at any height
- Travels virtually straight up and down
- Meets proposed entrapment guidelines
- Variety of assist / side rail options
- 500 lb. weight capacity

Benefits of Fast Rising

- Increases staff productivity
- Reduces staff injuries related to lifting
- Contributes to reducing staff turnover
- Saves staff energy and reduces accumulative trauma
- Meets OSHA speed recommendation
- Saves money




QD 8250



Transport / storage dollies

Options:


- Assist rail
- 1/2, 3/4, and full length side rails
- Side rail pads
- Orthopedic grid deck or perforated pan
- 360° mobility at all heights
- Main frame mounted trapeze
- Bed extender
- Individual motor lockout
- Foot end controls
- MX-Lock
- Variety of head and foot boards



MX-Lock

Specifications:	QD 8250	QD 8250Y
• Length (including bumper)	88"	84"
• Length (deck)	80"	76"
• Width (deck)	35"	35"
• Deck height - low	7 7/8"	7 7/8"
• Deck height - high	30"	30"
• Weight capacity	500 lbs.	500 lbs.

Specifications and dimensions are approximate and subject to change without notice • Patents pending



A wide variety of pressure reduction mattresses for MAXXUM series beds

Encore Medical Technologies, Inc.

12804 S Spring Street, Los Angeles, CA 90061
 (800) 558-8636, (310) 538-9138, (310) 538-9177 fax
www.encoremedical.com, email: info@encoremedical.com

