

SHOWER ACCESSORIES



Shower Stool S110

- Made of high impact plastic
- Compact size
- Large non-skid, suction rubber stops



Shower Bench S300

- Anodized aluminum frame
- Adjustable Seat Height
- Blow molded plastic seat with drainage holes
- 1" height increments



Bathtub Transfer Bench S650

- Anodized aluminum frame
- Padded vinyl seat with 1.5" foam
- Three piece seat
- 1" height increments



Transport Shower Chair S800

- Anodized aluminum frame
- Fits over most toilets
- Padded vinyl seat with commode pail
- Waterfall plastic armrests
- 4" swivel caster wheels with rear locks

SHOWER ACCESSORIES

Encore Medical offers an assortment of affordable brand shower accessories by Tuffcare. Please see individual specification information below.

SHOWER PRODUCT SPECIFICATIONS:

Shower Stool S110

3" High, 30" Wide, 16" Deep
Weight Limit: 260 lbs.

Shower Bench S300

14-19" High, 20" Wide, 12" Deep
Weight Limit: 260 lbs.

Bathtub Transfer Bench S650

19" – 23" High, 26" Wide, 16" Deep
Weight Limit: 250 lbs.

Note: Bariatric model available

Transport Shower Chair S800

20" High, 18" Wide, 16" Deep
Weight Limit: 250 lbs.



12804 S. Spring Street
Los Angeles, CA 90068
(800)_800-558-8636 Fax (310) 538-9177