

**M**illennium line

# Advanced Patient Care

MODEL ST370



FUNCTIONALITY

VERSATILITY

DURABILITY

ECONOMY



## Features - Benefits

- *the hydraulic shower trolley is a major cost saver enabling the washing and bathing of the most severely disabled patients to be carried out by only one patient nurse or attendant.*
- *twin-pedestal hydraulic height adjustment provides maximum stability in every situation and safe for even obese patients.*
- *height adjustment pedals are located on both sides of the bed.*
- *separate pedals serve for lowering the trolley - only light pressure needs to be applied to lower the platform gently and safely to comfortable low position - from either side.*
- *the mattress tray can be tilted to one side for draining and cleaning - the tilting mechanism is securely locked below the tray and is weight compensated by pneumatic gas spring - so requires no physical strength from attendants to operate.*
- *4 x 200 mm (8") castors, with central-locking brakes, and one steering castor for perfect mobility, with brake-pedals located at each corner.*
- *the shower tray is waterproof and gently padded, and is supplied with a removable head support pillow.*
- *the interior is manufactured from chemical cross-link 20 mm thickness PA 300 EVA covered with heat sealed PVC coated polyester - the foam repels water and cannot become waterlogged even in the event of the cover splitting.*
- *quick-release weight compensated drop-down side panels enable patients to be moved in and out of the trolley safely and quickly.*
- *the sturdy ABS undertrolley housing - to protect the mechanism from damage and rusting - includes built-in storage space for patients effects as well as an oxygen cylinder cradle, enabling the trolley to be used for bathing patients requiring oxygen assistance.*
- *height adjustment range 64 - 95 cm*

## Construction

*Twin-pedestal construction - the safest and surest way to ensure that the patient surface is always rigid and stable in every medical situation. The base frame manufactured from durable and robust galvanized steel tubes - is covered by a sturdy and elegant moulded ABS cover - which serves also as handy storage rack if necessary.*

*The mattress platform is constructed from HPL or fiberglass sheet reinforced with steel bars for absolute safety, and the entire chassis is finished in epoxy-polyester powder coating - in a wide choice of colour finishes.*

## Mobility

*The Millennium ST370 shower trolley is fitted with 8" (200 mm) castors, with 4-wheel locking and built-in steering mechanism which can be operated from any of the 4 corners of the stretcher by clearly indicated rocker pedals. The castors have thermoplastic rubber treads - ensuring longer life and smooth-running as well as improved mechanical behavior such as wear resistance, resistance to temperature extremes and improved mobility.*

## Colour Options

*The standard colour is blue (RAL 5010). other colours are available on request.*

# Millennium line

MODEL ST370

*Savion's Millennium Line new generation of patient care equipment - including hospital stretchers, ICU beds and shower trolleys - is designed for maximum patient comfort and safety - with minimum attendant effort and exertion*



## Standard features



*Central locking 8" castors brake/steer pedals at each corner.  
Drainage outlet and flexible drainage hose.*



*Fold-under side rail release lever.*



*Separate clearly marked pedals for raising and lowering located on both sides of the stretcher.  
Oxygen cylinder storage cradle.*



*Safety lock release lever for drainage and cleaning tilting function.*

# Specifications - Shower trolley model ST 370

## Overall

length:	202 cm	80"
overall width:	74 cm	29"

## Internal:

length:	180 cm	71"
width:	60 cm	23.5"

## Height Adjustment Range:

minimum height with 8" castors	64 cm	25"
maximum height with 8" castors	95 cm	37.5"

## Castors:

8" with central locking  
one tracking  
4 brake/steer pedals

## Options

- The shower tray can be placed in absolutely horizontal position for bathing treatment for burns cases by adjustment pin located below the mattress surface.
- All stainless steel castors.
- Stainless steel side rails.



\*Savion reserves the right to change specifications without notice



MODEL ST 360

# *Hydraulic Shower Trolley*



**SAVION**  
ADVANCED PATIENT CARE

## GENERAL

The new ST360 - another example of SAVION'S policy of constant improvement to its products - is an invaluable asset in cutting down labour costs in patient care. The pedal operated hydraulic height adjustment mechanism allows the trolley to be readily adjusted to the height of the patient. The quick-release drop-down sides provide full-length access for the patient to be placed in the trolley directly from the bed. The easy-running twin-wheel castors mean that the ST360 can be effortlessly taken into the shower room for patient bathing.

This whole operation from start to finish requires only one nurse or attendant - a major cost-saving factor in this type of patient care.

## FEATURES, OPTIONS AND ACCESSORIES

- sturdy all-steel chassis with reinforced fiberglass shower tray platform
- adjustable height by pedal operated hydraulic pump which can be operated from both sides of the trolley
- improved height adjustment range: 50 to 90cm
- waterproof padded shower tray and removable head support
- quick-release drop-down side panels

to enable patients to be moved in and out of the trolley quickly, easily and safely

- smooth-rolling 120mm twin-wheel castors in non-rusting plastic fitted with kick-stop brakes
- flexible drainage hose - the ST360 needs no plug or permanent drainage or installation



## CONSTRUCTION

The all-steel tubular chassis, and reinforced fiberglass shower tray platform ensure that the ST360 is strong and long-lasting. Castors - 120mm twin-wheel with all plastic construction - will not rust, and their easily adjustable kick-stop brake mechanism which locks both the swivel and the wheel itself, ensures that the trolley remains absolutely stationary whilst in use.

The side panels lift-up and drop-down to permit easy patient transfer, to and from the trolley - in complete safety and with minimum effort by the patient's attendant.

## HEIGHT ADJUSTMENT

The ST360 offers height adjustment variation of almost 40cm (15") - from a low position of 50cm (20") to a maximum of 90cm (35").

A further improvement is in the positioning and design of the height adjustment pedals, which provide a much narrower overall framework for the trolley, making it far more maneuverable in tight spaces and narrow passageways and doors. The overall width is now only 80cm (32") with the sides up, and less than 70cm (28") with both sides down.

## FINISH

Finished in epoxy paint, electrostatically applied and oven-baked, the ST360 has a rich and pleasant appearance. The epoxy finish is resistant to chipping and scratching. A choice of epoxy colours is available, with royal blue being standard for this product.



## TABLE OF DIMENSIONS

Specifications are subject to change from time to time

MATTRESS SUPPORT SURFACE				OVERALL			FWLER POSITIONING ANGLES	TRENDLENBURG AND REVERSE TRENDLENBURG ANGLES	
74"	22"	35"	20"	82"	closed 32"	open 28"	NONE	NONE	NONE
180 cm	57 cm	90 cm	50 cm	200 cm	80 cm	70 cm			

**Encore Medical Technologies, Inc.**  
 12804 S. Spring Street, Los Angeles, CA 90061  
 P: 800-558-8636, 310-538-9138, F: 310-538-9177  
[www.encoremedical.com](http://www.encoremedical.com),  
 Email: [info@encoremedical.com](mailto:info@encoremedical.com)