



**T.H.E. Medical**  
Professional Patient Care Products



## STEADY-AID 3500X

ELECTRIC  
SIT/STAND  
LIFT

- KEY FEATURES:**
- POWER BASE WIDENING
  - UP TO 450LB CAPACITY
  - REMOVABLE FOOTPLATE & KNEEPAD
  - OPTIONAL INTEGRATED WEIGH SCALE
  - 4-FUNCTION HAND CONTROLLER
  - LIFETIME FRAME WARRANTY

# STEADY-AID 3500X

## ELECTRIC SIT/STAND LIFT

THE INDUSTRY STANDARD FOR SAFETY & COMFORT



The Steady-Aid 3500X Electric Sit/Stand Lift is specifically designed to perform super-easy seat-to-seat transfers. It is ideally suited to assisting caregivers in the following tasks; peri-care, toileting, general transfers and rehabilitation of residents who are cognitive and co-operative. The Steady-Aid incorporates a power base widening system for optimal caregiver injury reduction. This feature eliminates arm and shoulder injuries caused by lifts with older style manual shifter levers. Additional features include a removable footplate and knee pad with calf strap, optional integrated digital weigh scale and 4 harness attach points for maximum versatility. An optional upgrade package increases the lifting capacity to an industry leading 450 lbs. An ergonomic control panel with a dual low battery warning system and dual controls means the Steady-Aid is always ready to use.



### Ergonomic Control Panel

- 4 function hand controller
- extra set of controls
- audio and visual low battery indicators
- emergency stop



### Removable Foot Plate and Knee Pad

- allows for rehabilitation and gait training



### Leg Support Harness

- provides additional support
- ideal for patients with minimal leg muscle tone



### Calf Strap

- prevents patients legs from slipping out of knee pad



# ADDITIONAL FEATURES AND OPTIONS



## REMOVABLE BATTERY PACK OPTION

Requires the use of a wall mounted charger and an extra battery pack.



## DIGITAL WEIGH SCALE OPTION

Allows for accurate patient weighing in lbs. and kgs.

## OUR COMPLETE FAMILY OF INTEGRATED LIFTING AND BATHING PRODUCTS



Sling Lifts



Bariatric Lifts



Ceiling Lifts



Recumbent Tubs



Side Entry Tubs



Reclining Tubs

## ORDERING INFORMATION

TSA3500X	Steady-Aid 3500X with On-Board Batteries and Charger
TSA3500XP	Steady-Aid 3500X with Removable Battery Pack and Wall-Mounted Charger
TSA3500XS	Steady-Aid 3500X with Integrated Digital Weigh Scale
TSA3500XPS	Steady-Aid 3500X with Removable Battery Pack and Integrated Digital Weigh Scale
T3410P	Removable Battery Pack - Extra
T10000	Low Base Caster Option (3" Front, 4" Rear-Braked)



**T.H.E. Medical**

Professional Patient Care Products

BARRIE, ON CANADA / LONGWOOD, FL USA

705-733-0022 FAX: 705-733-3432 CANADA/USA: 1-800-565-7075

WWW.THEMEDICAL.COM THESALES@THEMEDICAL.COM