



**T.H.E. Medical**  
Professional Patient Care Products



**ULTRALIFT  
3500X  
MOBILE SLING LIFT**

**KEY FEATURES:**


- 600LB LIFTING CAPACITY
- POWER BASE WIDENING
- EXTRA LONG BOOM
- SUPER 6 POINT SPREADER BAR
- 4-FUNCTION HAND CONTROLLER
- LIFETIME FRAME WARRANTY



# ULTRALIFT 3500X

## MOBILE SLING LIFT

THE INDUSTRY STANDARD FOR SAFETY & COMFORT



The Ultralift 3500X is the latest version of the industry's safest, most user-friendly mobile electric patient lifter ever. Manufactured in North America, the 3500X includes a long list of standard features specifically engineered to optimize patient comfort and safety, reducing caregiver injury. The 600 lb. lifting capacity allows for the safe transfer of your heaviest patients. The new power base widening system means the end of arm and shoulder injuries caused by old style manual shifter levers. The extra long boom eliminates patient/lift contact and allows for unobstructed 360 degree patient rotation. Floor lifts can now be accomplished "head first" without lifting the patient's legs over the main base or the base rails of the lift! Our signature "super six-point overhead spreader bar system" continues to provide for easy patient positioning while eliminating patient/spreader bar contact. Combined with our new repositioning slings, the Ultralift now performs patient rolling and in-bed repositioning at the push of a button.



### Ergonomic Control Panel

- 4 function hand controller
- extra set of controls
- audio and visual low battery indicators
- emergency stop



### Super 6-Point Spreader Bar

INCLUDES:

- pivoting, permanently attached side bars
- sling safety locks
- padded spreader bar



### Extra High Lifting Range

- combine the ultralift with the new walking harness for patient rehabilitation



### Extra Long Boom

- eliminates patient/lift contact
- significantly reduces injuries



### Super Easy "Head First" Floor Lifts

- eliminates all manual lifting



**T.H.E. Medical**  
Professional Patient Care Products

# ADDITIONAL FEATURES AND OPTIONS

## DIGITAL WEIGH SCALE OPTION

---



Allows for accurate patient weighing in lbs. and kgs.

## REMOVABLE BATTERY PACK OPTION

---



Requires the use of a wall mounted charger and an extra battery pack.

## WISHBONE STYLE SPREADER BAR

---



Not recommended for large and/or aggressive patients.  
Requires the use of wishbone style slings.

## OUR COMPLETE FAMILY OF INTEGRATED LIFTING AND BATHING PRODUCTS



Sit Stand Lifts



Bariatric Lifts



Ceiling Lifts



Recumbent Tubs



Side Entry Tubs



Reclining Tubs

## ORDERING INFORMATION

---

TUL 3500X	Ultralift 3500X with On-Board Batteries and Charger
TUL 3500XP	Ultralift 3500X with Removable Battery Pack and Wall-Mounted Charger
T3410P	Removable Battery Pack - Extra
TUL810	Optional Digital Weigh Scale
T10000	Low Base Caster Option (3" Front, 4" Rear-Braked)
TUL 2.5-310	Wishbone Style Spreader Option



**T.H.E. Medical**

Professional Patient Care Products

BARRIE, ON CANADA / LONGWOOD, FL USA

705-733-0022 FAX: 705-733-3432 CANADA/USA: 1-800-565-7075

WWW.THEMEDICAL.COM THESALES@THEMEDICAL.COM