

## STANDARD WHEELCHAIRS



### MODEL 227

- Folding powder coated frame with fixed full arms
- Embossed flame retardant vinyl upholstery
- Adjustable toggle wheel locks
- Cam Action swing away detachable footrests or leg rests
- Adjustable length foot plates
- High impact composite side panels

PLEASE SEE SPECIFICATIONS BELOW



### MODEL 277

- Folding triple chrome plated frame with detachable desk length arms
- Embossed flame retardant vinyl upholstery
- Adjustable toggle wheel locks
- Cam action swing away detachable footrests or leg rests
- Adjustable length foot plates
- Stainless side panels and high impact composite foot plates

PLEASE SEE SPECIFICATIONS BELOW



### MODEL 877

- Folding triple chrome plated frame with detachable desk length arms
- Padded double embossed flame retardant vinyl upholstery
- Padded arm and calf pads
- Adjustable toggle wheel locks
- Cam action swing away detachable footrests or leg rests
- Adjustable length foot plates
- Stainless side panels and aluminum foot plates

PLEASE SEE SPECIFICATIONS BELOW

## STANDARD WHEELCHAIRS

Encore Medical offers an assortment of affordable brand new wheelchairs by Tuffcare. Please see individual specification information below.

### WHEELCHAIR PRODUCT SPECIFICATIONS:

**Model 227:**

16" seat available with 250 lb weight capacity  
24" x 1" Rear molded wheels with hard rubber tires  
8" x 1" caster wheels with hard rubber tires  
Black upholstery

**Model 277:**

16" seat available with 250 lb weight capacity  
24" x 1" Rear molded wheels with hard rubber tires  
8" x 1" caster wheels with hard rubber tires  
Blue upholstery

**Model 877:**

16", 18" or 20" seat  
24" x 1" Rear molded wheels with hard rubber tires  
8" x 1" caster wheels with hard rubber tires  
Black upholstery  
Weight capacity 16" & 18" – 250lbs. 20" 310lbs.



12804 S. Spring Street  
Los Angeles, CA 90068  
(800)\_800-558-8636 Fax (310) 538-9177