

ELECTRIC BEDS

HILL-ROM 835/837 UPGRADE



Featuring:

Siderail Upgrade to Reduce Gaps

New White Finish

Patient Rail Controls

Tuck-away Side Rails

Central Brake & Steer System

(800) 558-8636 toll free
(310) 538-9138 local
(310) 538-9177 fax
www.encoremedical.com



ELECTRIC BEDS

HILL-ROM 835/837 UPGRADE

The Hill Rom 835/837 Upgrade is ideal for express recovery units. The functions include, head, knee and hi-low elevation along with cardiac chair and Trendelenberg positions. The four split side rails with upgrade inserts provide complete safety and can tuck away underneath the frame with a single push function. The 835/837 Upgrade can be controlled with dual side Central Brake & Steer Pedals and offers Automatic Contour at the knee section. This bed includes new plastic head & footboards and the Encore Medical 1 Year Warranty.

Includes: Central Break & Steer, Automatic Contour, Split Side Rails with Controls, Foot End Nurse Lockout Controls, 5" Casters, Spring or Pan Deck Surface and New Head & Footboards.

Options: IV Pole, Full Length Side Rails, ¾ Length Side Rails, Wall Bumpers, Trapeze Bracket, Pendant Hand Control, 3" or 4" Casters, Side Rail Communication.

HILL-ROM 835/837 UPGRADE DIMENSIONS:	HILL-ROM 835/837 UPGRADE FUNCTIONS:
Length – 88"	Head Up / Down
Width (Rails up) – 41.5"	Hi-Low Elevation
Width (Rails down) – 36.5"	Knee Up / Down
Highest – 30"	Cardiac Chair
Lowest – 18"	Trendelenberg
Sleep Surface – 35" x 80"	Reverse Trendelenberg
Weight Capacity - 500 lbs.	
Under bed Clearance – 7.5"	



(800) 558-8636 toll free
(310) 538-9138 local
(310) 538-9177 fax
www.encoremedical.com