

ELECTRIC BEDS

HILL-ROM CRITICAL CARE 894/896



Featuring:

Patient Rail Controls

Tuck-away Side rails

Central Brake & Steer System

CPR Release

Radiolucent Head & Foot Section

Patient Scale

Electric or Manual Operation

(800) 558-8636 toll free
(310) 538-9138 local
(310) 538-9177 fax
www.encoremedical.com



ELECTRIC BEDS

HILL-ROM CRITICAL CARE 894/896

The Hill Rom Critical Care bed is ideal for Intensive Care or Critical Care units. The functions include head, knee and hi-low elevation along with cardiac chair and hand controlled Trendelenberg positions. The four split side rails provide complete safety and can tuck away underneath the frame with a single push function. Radiolucent head section allows for in-bed X-Ray and this bed also includes a patient scale. The Critical Care bed also includes a CPR release, plastic head & footboards and the Encore Medical 1 Year Warranty.

Includes: Patient Scale, Radiolucent Head & Foot Section, Patient Nightlight, CPR Release, Central Brake & Steer, Split Side rails with Controls, Foot End Nurse Controls, 8" Casters, Pan Deck Surface and Plastic Head & Footboards.

Options: Air Mattress, IV Pole, Wall Bumpers, Pendant Hand Control, Side rail Communication.

HILL-ROM 835/837 CRITICAL CARE DIMENSIONS:	HILL-ROM 835/837 CRITICAL CARE FUNCTIONS:
Length – 91"	Head Up / Down
Width (Rails up) – 38.5"	Hi-Low Elevation
Width (Rails down) – 34"	Knee Up / Down
Highest – 40.5"	Cardiac Chair
Lowest – 22.5"	Trendelenberg
Sleep Surface – 33" x 80"	Reverse Trendelenberg
Weight Capacity - 500 lbs.	Electric & Manual Operation
Under bed Clearance – 7.5"	



(800) 558-8636 toll free
(310) 538-9138 local
(310) 538-9177 fax
www.encoremedical.com