

HILL ROM AFFINITY II

The Hill Rom Affinity Series birthing bed has been the industry workhorse in birthing units throughout the world. The functions include head, knee and hi-low elevations. This bed includes a refinished wood headboard and the Encore Medical 1 Year Warranty.

Includes: Placenta Basin, Calf Supports, Central Brake & Steer System, Nightlight, Lockout Controls, Instant CPR Release and New Custom fit Foam Mattress.

Options: IV Pole, Wall Bumpers, Side rail Communication, Pendant Hand Control, Air Mattress.

HILL ROM AFFINITY II DIMENSIONS:	HILL ROM AFFINITY II FUNCTIONS:
Length – 90”	Head Up / Down
Width (Rails up) – 42”	Hi-Low Elevation
Width (Rails down) – 34”	Knee Up / Down
Highest – 38”	Cardiac Chair
Lowest – 22”	Trendelenberg & Reverse Trendelenberg
Sleep Surface – 34 x 78”	
Weight Capacity - 500 lbs.	



12804 S. Spring Street
Los Angeles, CA 90068
(800) 800-558-8636 Fax (310) 538-9177